

Maintenance check of automated external defibrillator located on the wall of the Football Pavilion – 21st September 2020 – Clerk and Cllr Lynch

The following tests were carried out:

- ✓ Rescue Ready Indicator – Green, not action needed. The AED is ready to use (if the indicator, which can be viewed through the viewing window of the outer box, was red, refer to user guide).
- ✓ Check the voice prompts work and the display is readable
- ✓ Check that the defibrillation pads are ready for use and the service beep sounds
- ✓ Check that the LEDs work
- ✓ Check that the buttons work
- ✓ Check the AED case for signs of stress

A replacement battery was fitted to the alarm inside the external box and tested.

A new battery will need to be purchased by April next year for the AED. The cost of this would be approximately £300 (currently) inclusive of vat.

A monthly visual check could be carried out and recorded whilst the Play equipment inspection is carried out (to check the Rescue Ready Indicator is green not red. A six monthly check will be carried out by the Clerk and Cllr Lynch.